

GeneBlaster™
INSTRUCTION MANUAL

GeneBlaster



Gene



GeneBlaster Kits are new formulations of chemicals that significantly enhance and prolong the gene expression level obtained by *in vitro* transfection.

List of GeneBlaster Kits

Catalog Number	Description	Volume
GB20010	GeneBlaster Selection Kit ¹	3 x 1.5 mL
GB20011	GeneBlaster Ruby	3 x 1.5 mL
GB20012	GeneBlaster Sapphire	3 x 1.5 mL
GB20013	GeneBlaster Topaz	3 x 1.5 mL

¹ Contain 1 vial of each GeneBlaster reagent (Ruby, Sapphire & Topaz)

Use the content of the table above to determine the appropriate catalog number for your needs. You can order these products by contacting us. For all other supplementary information, do not hesitate to contact our dedicated technical support (tech@ozbiosciences.com).

OZ BIOSCIENCES

Parc Scientifique et Technologique de Luminy
Bâtiment grand luminy technopole
Zone entreprise – case 922
13288 Marseille Cedex 9
France

Tel: +33 (0)4.86.94.85.16

Fax: +33 (0)4.86.94.85.15

E-mail: contact@ozbiosciences.com

Web Site: www.ozbiosciences.com



OZ BIOSCIENCES
The art of delivery systems

1. Technology

1.1. Description

Congratulations on your purchase of the **GeneBlaster** Kit!

The **GeneBlaster** Kits are new formulations of chemical mixtures that significantly improve the gene expression level obtained with viral and non-viral gene delivery systems such as Magnetofection™, DreamFect™, DreamFect™ Gold, VeroFect and EcoTransfect reagents. Furthermore, they can be used with all commercially available transfection reagents. Since the application of the GeneBlaster Kits is cell type and promoter dependent three special formulations have been developed accordingly. They are extremely easy to use; simply add the appropriate GeneBlaster reagent to your culture medium and you get higher and longer levels of transgene expression.

- Rapid and easy to use
- Highest gene expression in many cells
- Prolong in vitro gene expression
- Excellent enhancement with transfection reagents such as Magnetofection™ and EcoTransfect™
- Effective for both transient and stable transfection
- Economical

1.2. Available Kits

OZ Biosciences offers four types of ready-to-use GeneBlaster Kits:

- 1. GeneBlaster Ruby** is a mixture of chemicals that has been developed intentionally for adherent cells. The dimension of the response is also cell type depend.
- 2. GeneBlaster Sapphire** is another combination of chemicals developed for adherent cells. The response extend is also cell type depend and complement well the GeneBlaster Ruby.
- 3. GeneBlaster Topaz** is a mixture of chemicals that has been developed purposely for suspension cells especially hematopoietic cell lines. However, other cell types appear to respond also very well.
- 4. GeneBlaster Selection Kit.** This kit is a convenient assortment of the three GeneBlaster reagents that permit to cover a large number of suspension and adherent cells.

1.3. Kit Contents

Each vial of the GeneBlaster Reagent (1.5 ml) is provided at a 100 X concentration and is sufficient for 150 transfections using 1 ml transfection volume. Each Kit contains three vials of reagents and allows performing at least 450 assays.

Kit contents

Description	Vials	Volume	Tube color
GeneBlaster Selection Kit ¹	3	3 x 1.5 mL ¹	1 red, 1 blue, 1 yellow
GeneBlaster Ruby	3	3 x 1.5 mL	3 red
GeneBlaster Sapphire	3	3 x 1.5 mL	3 blue
GeneBlaster Topaz	3	3 x 1.5 mL	3 yellow

¹ Contains 1 vial of each GeneBlaster reagent (Ruby, Sapphire & Topaz)

Stability and Storage

<u>Storage</u>	Upon receipt and for long-term use, store all reagent tubes at -20°C or + 4°C. GeneBlaster Kits are stable for at least 6 months at + 4°C or -20°C.
<u>Shipping condition</u>	The GeneBlaster Kits are shipped with gel pack (+4°C)

2. Applications

GeneBlaster Kits help to achieve higher and longer levels of transgene expression. The GeneBlaster Kits are new formulations of chemical mixtures that significantly improve the gene expression level obtained by transfection and infection. Since the GeneBlaster reagents are a mixture of specific chemicals; they could affect cell phenotype to a certain degree. These reagents might increase cellular toxicity for certain cells and optimization may be required (see §3.3 Important remarks).

Since the application of the GeneBlaster Kits is **cell type dependent** three special formulations have been developed. Indeed, some cells will respond greatly to one GeneBlaster formulation whereas others are not responsive. In another way, some cells might respond to the trans-activation induced by all or several chemical mixtures (see tables below).

In the same way, the effects of the GeneBlaster Kits are also **dependent on the nature of the promoter** that control gene expression (see table 3). Indeed, according to the promoter (CMV, SV40 early gene promoter, ubiquitin, EF1a, HIV LTR, or tissues/genes specific promoters, etc.) the effect might differ in particular cell lines.

2.1 Cell Types & Promoters

The gene expression enhancement and persistence mediated by the GeneBlaster reagent is cell type dependent and promoter dependent. These formulations have been successfully tested on a variety of suspension and adherent cells (see table 1 & 2). If a particular cell type or cell line is not listed below, this does not mean that GeneBlaster Kits is not going to work. An updated list of cells successfully tested is available on our website: www.ozbiosciences.com.

Likewise, these formulations have successfully trans-activated a variety of gene expression under the control of particular promoters (see table 3). If a particular promoter is not listed below, this does not mean that GeneBlaster Kits is not going to work. An updated list of responsive promoters is available on our website: www.ozbiosciences.com

Table 1: Cell Types with increased gene expression levels with the GeneBlaster Reagents

GeneBlaster Reagents successfully increase the expression level of a β -galactosidase, GFP, Luciferase or others genes in the cell lines listed. Promoters successfully tested were CMV, SV40 early gene promoter, EF1a, HIV LTR, or glucocerebrosidase-induced promoter (see Table 3).

Cell Line	Cell Type	Source	Ruby	Sapphire	Topaz
293, HEK-293, 293-T	Transformed Embryonic Kidney	Human	√	√	√
A549, H538, H460	Non-small cell lung carcinoma	Human	√	√	√
A172, F-98	Glioblastoma	Human	√	√	√
ALL	Acute lymphoblastic leukemia	Human	N/A	N/A	√
Caco-2	Colon Adenocarcinoma	Human	√	√	N/D
CHP126	Neuroblastoma	Human	√	√	N/D
CHO	Epithelial-like (Ovary)	Hamster	√	√	N/A
HeLa, Hela-S3	Cervical Epithelial Carcinoma	Human	√	√	N/A
HT 29	Colon Adenocarcinoma	Human	√	√	N/D
HUVEC	Endothelial Cells (primary)	Human	N/A	√	N/A
Jurkat	T leukemia cells	Human	N/A	N/A	√
K-562	Myelogenous Leukemia	Human	N/A	N/A	√
MCF-7, Hs578T	Breast Adenocarcinoma	Human	√	√	N/A
NIH3T3	Fibroblasts	Mouse	√	√	N/A
P19	Teratocarcinoma	Mouse	√	√	√
PC-12	Pheochromocytoma (adrenal)	Rat	√	√	N/A
Primary	Peripheral Blood Lymphocytes Peripheral Blood Mononuclear cells	Human Mouse	N/A	N/A	√
Primary	Neurons	Rat	√	√	N/D
S2 TnB1-4	Drosophila melanogaster Trichoplusia ni BTI	Insect	√	√	N/D

LEGEND: √: Works, N/A: Not Appropriate (due to toxicity and lowered gene expression) and N/D: Not Determined

Table 2: Gene expression controlled by CMV promoter in Cells Non-Responding to the GeneBlaster Reagents

For some of the cell lines listed, GeneBlaster may exhibit some cytotoxicity. In these cells, only a CMV promoter has been analyzed. Since GeneBlaster effect is promoter dependent, other promoter might be trans-activated by the chemical cocktails.

Cell Line	Cell Type	Source
B16F10	Melanoma	Mouse
CHO-K1	Epithelial-like (Ovary)	Hamster
COS-7	Fibroblast (Kidney)	Green Monkey
CV-1	Fibroblast-like (Kidney)	Monkey
HepG2	Hepatoma	Human
MDCK	Normal -Kidney	Canine

Table 3: Promoters sensitive to trans – activation mediated by the GeneBlaster Reagents

The effects of the GeneBlaster Kits are also dependent on the nature of the promoter that controls gene expression. The nature and degree of trans activation depend also the cell types used (see tables above).

Promoters	Ruby	Sapphire	Topaz
CMV	√	√	√
SV-40 Early gene promoter	√	√	√
EF1a	√	√	√
HIV LTR, MoMSV LTR, RSV LTR...	√	√	√
MAT2A, IL2, MIP1alpha, IFN γ , NF-kappa B...	N/D	N/D	√
M human ChAT, H1(0), CTMP, MMTV, GPH alpha, Cytochrome P450, PLAP-1, Glucocerebrosidase-induced ...	√	√	N/D

LEGEND: √: Works, N/A: Not Appropriate (due to toxicity and lowered gene expression) and N/D: Not Determined

3. Protocols

3.1. General Considerations

- Select the suitable GeneBlaster reagent for favored cell types according to the Table 1.
- Briefly vortex the reagent before each use. If stored at -20° C, bring up to room temperature.
- Dilute the GeneBlaster reagent 100 times in the culture medium.

The instructions given below represent sample protocols that were applied successfully with a variety of cell lines and promoters. They can be used as guidelines to achieve very high gene expression level with minimal times. Optimal conditions do vary from cell line to cell line, promoter to promoter and the final dilution of the GeneBlaster Reagents might have to be adjusted to achieve best results. Therefore, we advise you to optimize few trans-activation parameters (concentration, incubation time, medium change...) if necessary.

3.2. General Protocol

1. Prepare the DNA/transfection reagent complexes (such as Magnetofection™ DreamFect™ or EcoTransfect reagents) according to the manufacturer's instruction or the viral titers as standard.
2. Add the complexes or viruses onto the cells growing in serum-free or serum-containing medium as standard culture conditions or as suggested by the manufacturer's protocol.
3. Incubate 4 hours.
4. Then, add the proper GeneBlaster to the cells (see tables 1, 3 & 4). Each GeneBlaster reagent supplied is 100X concentrated. The best method, to add the GeneBlaster to the cells in order to achieve the best mixing and to minimize cytotoxicity, is to first dilute the GeneBlaster in serum-

containing medium. Therefore, prepare serum-containing medium with GeneBlaster at 2 to 10 X concentration (i.e., 1:50 to 1:10 dilution). For example, if the cells are transfected in 1ml of medium, prepare 1 ml of 2X (1:50 dilution) or 250 µl of 5X (1:20 dilution) or 110 µl of 10X (1:10 dilution) GeneBlaster solution in serum-containing medium. Finally, transfer the culture medium containing GeneBlaster to the cells (final concentration is 1X).

- **OPTION.** As an alternative, the GeneBlaster can be directly added to the cells at a 1:100 dilution (1X final concentration, i.e. 10 µl in 1 ml) and mixed immediately by swirling or very gentle pipetting. The mixing process is very important to avoid localized toxic effects. If the transfection is realized in serum-free medium, first add the serum-containing medium to the cells and next, add the GeneBlaster to 1X final concentration.
5. The next day and if required, add fresh growth culture media without GeneBlaster. If some toxicity is observed and to avoid potential toxic effects, 24 hours post-transfection proceed to a medium change (discard the GeneBlaster -containing medium and replace by a fresh complete growth medium without GeneBlaster.
 6. Cultivate the cells under standard conditions until evaluation of gene expression. The gene expression analysis can be monitored and assayed 24 to 72 hours following transfection or infection. This depends on the cell type, reporter gene and promoter activity.

Stable Transfection: the GeneBlaster reagents can also be used to produce stably transfected or transduced cells. 48 hours post-transfection or 24 hours post-infection replace old growth medium by fresh medium containing the appropriate antibiotics for selection.

Table 4: Recommended Uses for the GeneBlaster Reagents

Cells	Ruby	Sapphire	Topaz	Fold Improvement
293, HEK-293	++	++	++	2-3X
A549	++	++	N/D	4-5X
Caco-2	++	N/D	N/D	2-4X
HeLa	+++	++	N/A	5-10X
HT 29	+	N/D	N/D	4-5X
HUVEC	N/A	+	N/A	2-3X
Jurkat	N/A	N/A	+++	5-10X
K-562	N/A	N/A	+++	5-15X
MCF-7	++	++	N/A	2-10X
NIH3T3	++	+	N/A	N/D
P19	+++	+	++	2-3X
PC-12	++	+++	N/A	3-10X

LEGEND: + Works Well, ++ Works Better, +++ Works Best, N/A Not Appropriate (due to toxicity and lowered gene expression) and N/D: Not determined.

3.3. Important Remarks

- **For all GeneBlaster Kits**

Although the team of OZ Biosciences has carefully designed and optimized the GeneBlaster formulations for a number of cells, additional adjustment and optimization might be required for other cells to minimize toxicity and to enhance gene expression level after transfection or infection. We advise testing various dilutions of the GeneBlaster from 50X to 250X dilutions to attain the best results.

- **Caution!**

For some cell types such as B16-F0, BHK-21, CHO-K1, COS-7, CV-1, Hep-G2 & MDCK (especially if a CMV promoter controls the gene of interest) we do not recommend the use of the GeneBlaster Reagents. The chemicals either have no effect on the expression level of the gene controlled by a CMV promoter or on the transfection/infection efficiency or increase toxicity. Consequently, in function of the cell, the promoter and the concentration of GeneBlaster the chemicals mixture may induce inferior transfection efficiency and elevate toxic effects.

- **For GeneBlaster Ruby & Sapphire**

- ✓ For handiness, GeneBlaster Ruby & Sapphire can be transfer to the cell culture medium at different times during the procedure. For example, they can be added to the cells concurrently with the DNA/transfection reagent complexes (or virus) or 4 hours post-transfection/post-infection (see section 3.2). Contrary to Topaz, Ruby & Sapphire formulations are not efficient if added 24 hours after transfection.
- ✓ These two formulations can be relatively toxic for some cells. Consequently, the culture medium may need to be changed the next day with fresh GeneBlaster -free culture medium.

- **For GeneBlaster Topaz**

For convenience, GeneBlaster Topaz Reagent can be added to the cell culture medium at different times during the protocol. It can be added to the cells simultaneously with the DNA/transfection reagent complexes (or virus) or 4 hours post-transfection/post-infection (see section 3.2) or 24 hours after transfection/infection.

4. Related Products

Description	Reference
Magnetofection Technology	
Super Magnetic Plate	MF10000
Magnetic Plate 96-magnets	MF10096
PolyMag 1mL (<i>for all nucleic acids</i>)	PN31000
CombiMag 1mL (<i>to boost transfection reagent</i>)	CM21000
SilenceMag 1mL (<i>for siRNA application</i>)	SM11000
ViroMag 1mL (<i>for viral application</i>)	VM41000
ViroMag R/L 1mL (<i>for lenti- & retro-viruses</i>)	RL41000
NeuroMag 1mL (<i>for transfection of neurons</i>)	NM51000
SelfMag Amino Kit	SA10000
SelfMag Carboxy Kit	SC20000
Protein Delivery Systems	
Ab-DeliverIN 1 mL	AI21000
Pro-DeliverIN 1 mL	PI11000
Tee-Technology (lipid-based reagents)	
Lullaby siRNA transfection reagent	LL71000
DreamFect Gold Transfection reagent 1mL	DG81000
EcoTransfect Transfection Reagent 1mL	ET11000
VeroFect Transfection Reagent 1mL	VF61000
FlyFectin Transfection Reagent 1mL	FF51000
Gene & Protein Tools	
Bradford – Protein Assay Kit	BA00100
β-Galactosidase (ONPG) assay kits	GO10001
β-Galactosidase (CPRG) assay kits	GC10002
X-Gal Staining Kit	GX10003

Our dedicated and specialized technical support group will be pleased to answer any of your request and to assist you in your experiments. Do not hesitate to contact us for all complementary information and remember to visit our website in order to stay inform on our last breakthrough technologies and updated on our complete product list. <http://www.ozbiosciences.com>

Limited License

The purchase of the **GeneBlaster** Kit grants the purchaser a non-transferable, non-exclusive license to use the kit and/or its separate and included components (as listed in section 1, Kit Contents). This reagent is intended **for internal research only** by the buyer. Such use is limited to the use described in the product manual. In addition, research only use means that this kit and all of its contents are excluded, without limitation, from resale, repackaging, or use for the making or selling of any commercial product or service without the written approval of OZ Biosciences.

Separate licenses are available from OZ Biosciences for the express purpose of non-research use or applications of the **GeneBlaster** Kit. To inquire about such licenses, or to obtain authorization to transfer or use the enclosed material, contact the Director of Business Development at OZ Biosciences.

Purchasers may end this License at any time by returning all **GeneBlaster** Kit material and documentation to OZ Biosciences, or by destroying all GeneBlaster Kit components. Purchasers are advised to contact OZ Biosciences with the notification that a GeneBlaster Kit is being returned in order to obtain a refund and/or to expressly terminate a research only license granted through the purchase of the kit(s).

This document covers entirely the terms of the **GeneBlaster** Kit only license, and does not grant any other express or implied license. The laws of the French Government shall govern the interpretation and enforcement of the terms of this License.

Product Use Limitations

The **GeneBlaster** Kit and all of its components are developed, designed, intended, and sold for research use only. They are not to be used for human diagnostic or included/used in any drug intended for human use. All care and attention should be exercised in the handling of the kit components by following proper research laboratory practices.

For more information, or for any comments on the terms and conditions of this License, please contact:

Director of Business Development
OZ Biosciences
Parc Scientifique et Technologique de Luminy
Bâtiment grand luminy technopole – zone entreprise case 922
13288 Marseille Cedex 9
France
Tel: +33 (0)4.86.94.85.16
Fax: +33 (0)4.86.94.85.15
E-mail: contact@ozbiosciences.com