

EcoTransfect™: List of successfully transfected cells

Cell Lines			
Cell	Type	Species	% transfected cells
293	Embryonic kidney	Human	75
293T	Embryonic kidney	Human	75
3T6	Embryonic fibroblast	Mouse	40
A549	Non-small cell lung carcinoma	Human	35
B16F10	Melanoma	Mouse	45
BEAS-2B	Bronchial epithelium	Human	55
BHK21	Kidney	Syrian Hamster	OK
CHO	Ovary epithelial like	Chinese Hamster	85
CHO-K1	Ovary epithelial like	Chinese Hamster	85
COS1	Kidney	Green Monkey	80
COS7	Kidney	Green Monkey	80
CV-1	Fibroblast-like Kidney	Green Monkey	OK
FRT	Fisher Rat Thyroid	Rat	30
HEK293	Embryonic kidney	Human	75
HeLa	Cervix carcinoma	Human	50
MDCK	Kidney	Canine	35
N2A	Neuroblastoma	Mouse	55
NIH3T3	Fibroblast	Mouse	50
U87	Glioblastoma-astrocytoma	Human	50
Vero	Kidney	Green Monkey	35