



## Lullaby® - List of cells successfully tested

**OZ BIOSCIENCES**  
The art of delivery systems

Cell line	Cell type & description	
1207	Human carcinoma	+
293, 293T	Human epithelial	+
3T6	Mouse fibroblast	+
A-293	Human embryonic	+
A549	Human carcinoma	++
B16-F10	Mouse melanoma	+
B16-OVA	Mouse Melanona	+
BEAS-2B	Human epithelial	+
BHK-21	Human fibroblast	+
CHO, CHO-K1	Human epithelial	+
COS-1, COS-7	Monkey fibroblast	+
CV-1	Monkey fibroblast	+
H441	Human carcinoma	+
HEK293	Human embryonic	+
HeLa	Human adenocarcinoma	+
IC21	Mouse Macrophage	+
K562	Human Myelogenous leukemia	++
M-1	Mouse epithelial	+
MCF-7	Human adenocarcinoma	+
MDAMB435	Human Breast Cancer	+
MDCK	Canine epithelial	+
MEF	Embryonic fibroblast	+
MIA PaCa-2	Human carcinoma	+
MIN6	Mouse pancreatic insulinoma	+
MRC5	Human Lung Fibroblasts	++
N2A	Mouse neuroblastoma	+
NCI-H441	Human Lung adenocarcinoma	+
NIH-3T3	Mouse fibroblast	+
PC-12	Rat Pheochromocytoma	+
Pre-B	Mouse lymphoma	+
SH-SY5Y	Human neuroblastoma	++
U2-OS	Human Osteosarcoma	++
U251	Human Glioblastoma	++
U87	Human glioma	+
Vero E6	Monkey epithelial	++

Primary cell	Origine	
Dermal Fibroblasts	Human	++
Embryonic stem cells (hESC)	Human	++

+	successfully tested
++	successfully tested and published

*Lullaby® List of cells - March 2010*